Versuch Immunfluoreszenz

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Immunfärbung

Needed

Pasteur pipettes, Pasteur pipette bulb, forceps, Bunsen-burner, cover slips and microscope slides, Petri dishes, sterile Minimal medium.

Solutions

Fixative (92.6 ml)

200mM	PIPES	PH6.7 with 10N NaOH	50ml (prepare under exhaust hood!)
500mM	EGTA	PH8.5 with 10N NaOH	10ml (prepare under exhaust hood!)
1 M	MgSO ₄		1 ml
	DMSO		10 ml
	Formaldehyde (37%)		21.6 ml

10X Phosphate buffered saline (PBS) (1L),

80 g NaCl, 2 g KCl, 14.4 g Na₂HPO₄, 2.4 g KH₂PO₄
 With 800ml H₂O, PH 7.4, then fill up to 1000ml, autoclave.
 Prepare (1L) **1X PBS**

NaCitrat (50mM, pH 5.8), prepare 50 ml

NaCitrat (100mM, pH 4.5), prepare 50 ml

Digestion solution (5 ml)

Driselase	100 mg	(solved in 600µl 50 mM sodium citrate pH 5.8)
Zymolyase	20 µl	(5U/µl)
GlucanX	800 mg	(solved in 1 ml 50 mM sodium citrate pH 5.8)
β-D Glucanase	40 mg	(solved in 400 µl 100 mM sodium citrate pH 4.5), heating at 55°C for 5min, then incubating on ice 30min to reduce protease activity).

Collect all enzymes together and fill with 50 mM sodium citrate pH 5.8 to 2,5 ml. Add 2,5 ml egg white and vortex (**egg white** could be aliquoted and frozen at -20°C).

10X TBS (1L),

20mM Tris/HCl, pH 7.6	(24.2 g/l)
137mM NaCl	(80 g/l)

1X TBST (1L)

1XTBS + 0.1% Tween 20 (1 ml/l)

Blocking solution (100 ml)

1XTBST + 5% skim milk (5 g/100 ml)

Mounting media with DAPI from VECTASHIELD.

3l ddH₂O autoklavieren

Minimal medium for *A. nidulans* growth

50 ml Salt stock solution; 1 ml Microelement stock solution; 20 g Glucose; adjust to

pH 6.5 using 10 N NaOH (plus the auxotroph markers).

Antibodies

Antibody dilutions: First antibody: anti alpha tubulin (DM1A + TUB-1A2, 1:1). Second antibody: anti mouse IgG-FITC (Sigma).

Toxicology

Formaldehyde is irritant. When the stock solution is made, weighing, boiling and aliquotation have to be done under the fume hood. Wear glasses and gloves. It can etch your eyes and skin!

Methanol is extremely toxic. If ingested, as little as 10 ml can cause permanent blindness and as little as 60 ml can result in death.

Procedure for germlings

- place flamed and sterile cover slips on bottom of petri dish
- add 450 µl microscope medium (with glucose as C-source)
- Inoculate with 10^3 spore/ml, let it grow for 12-24 hours at RT (prepare therefore a 10µl spore suspension, add one drop of the suspension on each coverslip).
- Remove media and add **immediately** fixation solution, incubate at RT for **30 min** (not longer!)
- Three 10min washes in 1XPBS. In the last washing step the coverslip should be carefully dipped into the buffer, hanging buffer dripped into towels. Switch to a new petri dish.
- Incubate with 200 µl digestion solution for 1h at RT (Driselase is insoluble that's why you have to mix the digestion solution permanently before using. If necessary cover with parafilm on top of the cover slips).
- Three 10 min washes in 1XPBS at RT. In the last washing step the coverslip should be carefully dipped into the buffer, hanging buffer dripped into towels. Switch to a new petri dish.
- 10 min (**exactly**) incubation at -20°C in -20°C **precooled** methanol.
- Two 5 min washes in 1XPBS. In the last washing step the coverslip should be carefully dipped into the buffer. Switch to a new petri dish.
- 15 min incubation in blocking solution (1XTBST + 5% skim milk). Switch to a new petri dish.
- First antibody over night at 4°C , anti-tubulin 1:500 in (1XTBST + 5% skim milk), cover with parafilm on top of the cover slips.
- Four 10 min washes with 1XTBST. In the last washing step the coverslip should be carefully dipped into the buffer, hanging buffer dripped into towels.

Switch to a new petri dish.

- Second antibody at RT for 1 h. FITC antimouse 1:100 in TBST.
- Four 10 min washes with 1XTBST. In the last washing step the coverslip should be carefully dipped into the buffer, hanging buffer dripped into towels.
- Mount on microscope slide above a drop of mounting media (about 10 μ l).
- Seal with nail polish, it has to be completely sealed up! store at 4°C until viewing.

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Protoplastentransformation
by Fischer & Liese, adapted by Requena

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Minimal medium (MM)

50 ml Salt stock solution; 1 ml Microelement stock solution; 20 g Glucose; (for Protoplast add 0.6M KCl) adjust to pH 6.5 using 10 N NaOH.

20 x Salt stock solution

120 g NaNO₃; 10.4 g KCl; 10.4 g MgSO₄ x 7H₂O; 30.4 g KH₂PO₄.

1000 x Microelement stock solution:

22 g ZnSO₄ x 7H₂O; 11 g H₃BO₃; 5 g MnCl₂ x 4H₂O; 5 g FeSO₄ x 7H₂O; 1.6 g CoCl₂ x 5H₂O; 1.6 g CuSO₄ x 5H₂O; 1.1 g (NH₄)₆Mo₇O₂₄ x 4H₂O; 50 g Na₄EDTA; adjust to pH 6.5-6.8 using KOH

Table 3.6 Vitamins, amino acids and medium components

Component	Stock Concentration	Volume per liter
Pyridoxin-hydrochloride	0.1 %	1 ml
Uracil		1 g
Uridine		1.2 g

Microscopy Minimal medium (MM)

50 ml Salt stock solution; 1 ml Microelement stock solution; 20 ml Glycerol; adjust to pH 6.5 using 10 N NaOH.

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Everything under sterile conditions.

Clean the bench with 70 % ethanol, as well as your hands.

To open the flask: flame the neck and top of the closed flask, then open the lid and flame once more the opening of the flask.

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- Pre-culture from plate or conidiophore suspension (strain TN02A3)
- Inoculate liquid medium (500 ml) and incubate O/N (12-16 h) 30°C in the waterbath
- Filter the culture through sterile miracloth using the porcelain funnel. Wash with wash solution. Collect with a little spoon (steril!)
- Set 1g mycelium in a 100 ml sterile Erlenmeyer flask with 5 ml of osmotic medium (careful not to overexces the volume because then it is difficult that protoplasts can migrate through upwards)
- Add GlucanX, 180 mg dissolved in 1 ml ddH₂O and incubate on ice for 5 minutes.
- Add BSA (10 mg dissolved in 0.5 ml of sterile water), swirl and incubate in shacker for 1 (-2)h @30°C with gentle shaking 70 rpm. Check at the microscope
- Transfer to a sterile 50 ml falcon and overlay with 10 ml trapping buffer carefully and slowly close to the surface of the liquid
- Centrifuge 15 minute @ 5000 rpm, 4°C in HP-6 rotor (free swing), (be careful: slow accelerations!!)
- If you see a band of protoplast collect it with a sterile pasteur pipette and transfer to a sterile 15 ml corex tube and add STC (fill the whole tube up to 12 ml)
- Centrifuge @ 7000 rpm, 4°C for 10 minutes
- Remove the supernatant, add gently 2 ml of STC without resuspending the pellet and remove it again
- Resuspend the pellet in the remaining drop (ca. 100-200 µl). Take 10 µl and count the cells in the counting camera
- Dilute the protoplast to an appropriate concentration of 10⁶ prot/100 µl

Transformation

- 100 μ l pf protoplast (10^6)
- 100 μ l of DNA (approx. 10 μ g of DNA filled with 1x STC*)
* If vol of DNA > 10 μ l add same volume of 2x STC and fill up to 100 μ l with 1xSTC
- 25 minutes incubation at RT (in a falcon tube)
- Add 2 ml of PEG and roll the tube until it looks homogenous
- 20 minutes incubation @ RT
- Add 8 ml of STC and roll it again until total homogenization
- Spread different volumes on appropriate plates (MM+KCl) + appropriate markers (in our case no marker is needed), Incubate 3 days @37°C and then pick for the microscopy.

Microscopy

- Flamed and sterile cover slips on bottom of petri dish, add 0.4 ml MM glycerol on top
- With sterile toothpicks pick the colonies, and put some spores into media, another copy will be taken on small minimal media plate (after making marks with raster plate) grow for 12-24 hours at RT before microscopy

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